

# Test Protocol Outline

## Corrositex®: In-Vitro Corrosivity Test for Transport - Class 8

Determines chemical corrosivity and permits assignment of Packing Group classification for Class 8 corrosives. This test replaces the rabbit test of dermal corrosivity by providing a reliable means of mimicking this test. U.N. Packing Group classification can be determined or marketing claims substantiated.

### Experimental Design

The proprietary core technology of the **Corrositex** test is based on a biomembrane and chemical detection system.

A proprietary bio-barrier membrane is penetrated by the sample, thus mimicking the effect of corrosives on living skin. As the corrosive sample destroys this bio-barrier, the fluid below changes color or texture. This time for these changes to occur is determined and related to previously determined pH behaviour.

### Method

#### Step 1: Qualification

The sample is first screened to determine suitability for the test.

#### Step 2: Categorisation

Small aliquots of the test sample are placed into the **Corrositex** categorising testing system. The colour change gives the category assignment - Category I or II.

#### Step 3: Classification

Small aliquots of the test sample are placed into the **Corrositex** biobarrier system. This consists of a glass vial filled with a chemical detection fluid capped by a freshly prepared proprietary bio-barrier membrane. The membrane is designed to mimic the effect of corrosives on living skin. As soon as the corrosive sample destroys this bio-barrier, the fluid below changes color or texture. The time taken for penetration through the membrane is determined.

**Step 4: Assignment** The appropriate U.N. Packing Group classification can now be assigned...

| Category        | Time Required for CDS changes (minutes) |              |                |               |
|-----------------|---|--------------|----------------|---------------|
| Category 1      | 0 to 3 min                              | >3 to 60 min | >60 to 240 min | >240 min      |
| Category 2      | 0 to 3 min                              | >3 to 30 min | >30 to 60 min  | >60 min       |
| Packaging Group | I                                       | II           | III            | Non-corrosive |

### Recognitions

Department of Transportation - DOT-E 10904  
 Consumer Product Safety Commission (CPSC)  
 European Centre for the Validation of Alternative Methods (ECVAM)  
 EPA Federal Register / Vol. 60, No. 142 Dermal Corrosion  
 Food & Drug Administration (FDA)  
 International Air Transportation Association (IATA)  
 National Institute of Environmental Health Sciences (NIEHS)  
 Occupational Safety and Health Administration (OSHA)  
 Transport Canada - Permit For Equivalent Level of Safety SU 4483

A full list of specific references is available on request.

### References

In-vitro Corrositex Reference Manual DOT - E 10904